


# Calibrating pH Meters

## Get to Know Your Aquaseacher AB23PH OHAUS Bench Meter

Range: 0.00–14.00 pH  
Resolution: 0.01 pH  
Accuracy: 0.05 pH


Scan QR code or go to  
[knowledge.carolina.com/labskills](http://knowledge.carolina.com/labskills)  
for videos and more!



- Calibration button
- pH electrode
- Power button
- Temperature reading
- LED display
- pH reading
- Read/Enter button

# Get to Know Your OHAUS ST20 Pen Meter

**Range:** 0.00–14.00 pH  
**Resolution:** 0.01 pH  
**Accuracy:** 0.05 pH



- Calibration button
- Hold button
- pH electrode
- Power button
- Electrode protection cap
- Display
- pH reading
- Temperature reading